QUINTIP®

Individual Skin Test System

FOR PUNCTURE SKIN TESTING

QUINTIP® is a skin test device designed to apply test extract using the puncture technique.

Figure 1

QUINTIP® is designed for use by allergy practitioners who are trained in the application and interpretation of the puncture technique and who are trained in the recognition and treatment of adverse allergic reactions should they occur.

Figure 2

QUINTIP® System consists of three components: The **QUINTIP®** skin test device (Figure 1), the **ComforTen®** 60-hole covered tray (Figure 2), and the reservoir (Figure 3). A 30-hole tray is also available.


Figure 3

PRODUCT DESCRIPTION

Each **QUINTIP®** device has a molded plastic grip with one test probe (Figure 1). Inserted into each unit is a stainless steel lancet tip that protrudes from the probe enough to give the proper testing depth. The **QUINTIP®** device is to be used once and discarded in an approved sharps container. The sterile test devices are provided as 20 units per pouch with 15 pouches in a shelf pack for a total of 300 units.

The stackable, covered trays hold 30 or 60 reservoirs. They are designed for testing and storage flexibility. Both trays are provided as single units and include optional-use spacers (Figure 5) for increasing the number of tests between filling the reservoirs.

The reservoirs are supplied in sterile pouches of 30. The **QUINTIP®** device covers the reservoir, helping protect the antigen from foreign material and evaporation.

SET-UP INSTRUCTIONS

Remove the **ComforTen®** covered tray from the box. Each **ComforTen®** covered tray is supplied with blank self-adhesive labels and spacers (Figure 4). Each label strip has five designated areas corresponding to the reservoirs. Note the name of each allergenic extract in the area provided. Remove the label backing and place on the tray inline with the corresponding reservoirs.


Figure 4

HollisterStier Skin Test Reservoirs are provided in a sterile pouch of 30. Peel back the paper label as indicated on the pouch of the reservoirs. Firmly place a reservoir into each hole in the

ComforTen® tray. If using the optional spacers, fill the appropriate reservoir with approximately 20 drops (not more than 0.8 mL) using the dropper from the HollisterStier Allergy 5mL scratch test extract vial or from the bulk vial using a sterile syringe (Figure 5b). Place the spacers on top of the reservoirs as noted in Figure 5a. If the spacers are not used, fill the appropriate reservoir with approximately 12 drops (not more than 0.5 mL) using the dropper from the HollisterStier Allergy 5mL scratch test extract vial (Figure 5c). Store the unused spacers in the slots provided on the tray.


Figure 5a

Once the trays are prepared, remove a **QUINTIP®** device from the sterile pouch by peeling off the label as indicated on the package. Pick up the device by the plastic grip using caution not to touch any part of the leg or the lancet tip. Place the tips in the prefilled reservoirs. Continue until all of the reservoirs are filled.

leg or the lancet tip. Place the tips in the prefilled reservoirs. Continue until all of the reservoirs are filled.


Figure 5b

Figure 5c

The 30-hole and 60-hole trays hold 30 or 60 **QUINTIP®** devices respectively. With the **QUINTIP®** devices in the tray, the reservoirs are covered to prevent foreign matter from entering the extract. Once the tray is loaded, it can be used immediately or covered and stored at 2°- 8°C until needed.

SKIN TEST PROCEDURE

Prior to testing, clean the test area with alcohol and allow to dry. Refer to "Highlights of Prescribing Information" and "Positive Skin Test Control-Histamine" or other manufacturers' package inserts for more detailed testing instructions. Remove a **QUINTIP®** skin test device from a prepared tray by lifting it vertically, being careful not to touch the tip against the reservoir. With the **QUINTIP®** perpendicular to the skin, press down on the skin with medium pressure without lifting the device from the skin. Remove by lifting vertically and discard the **QUINTIP®** in an approved sharps container. Small circles about 4mm in diameter should remain at the test site indicating the correct amount of pressure was applied. The visibility of the circles will vary between patients according to the thickness, fragility, and pigmentation of their skin. The amount of pressure needed to produce a satisfactory test site requires that the depth guard come in full contact with the skin.

Following the above instructions will produce uniform skin tests all at the proper depth without inducing bleeding.

STORAGE AND DISPOSAL

Once a patient's test session has been completed, place new **QUINTIP®** devices in the antigen trays and store at 2°- 8°C until the next test session.

QUINTIP® devices are designed to be used once and discarded (in accordance with all local, state and federal regulations). These devices cannot be cleaned sufficiently to prevent cross contamination from repeated use or to protect subsequent patients from possible transfer of serum hepatitis or other bloodborne pathogens.

QUINTIP® reservoirs should be replaced periodically. For best results, discard and replace the reservoirs as they are emptied and refill with fresh extract. If a used **QUINTIP®** device is inadvertently placed in a reservoir, the individual reservoir should be discarded and a new reservoir prepared before the next test session.

READING SKIN TEST REACTIONS

Measure wheal and erythema responses for histamine positive control at 10 minutes and allergens at 15 and 20 minutes.² For optimal response, use a positive control with a 6mg/mL histamine base. Immediate reactions to histamine typically peak at 8 minutes while allergens peak at 15 minutes.¹ HollisterStier Allergy recommends recording wheal and erythema measurements to improve accuracy, precision and uniformity of diagnostic testing and eliminate confusion as to the quantitation of the allergic response.

Reactions from the **QUINTIP®** device may be graded as follows:

European Prick Test Grading System

Grade	Wheal (mm)	Erythema (mm)
Negative	<3	Not Measured
Positive	>3	Not Measured

Graduated System

Grade	Wheal (mm)	Erythema (mm)
0+	<3	<15
2+	<6	<15
3+	7-9	16-30
4+	>9	>30

A positive skin reaction to any allergenic extract must be interpreted in light of the patient's history of symptoms, known allergen exposures, and time of the year.

The skin tests are in no way a substitute for a careful allergic history. Rather, they serve as additional information to aid in identifying causative allergens in patients with allergic disorders.

HOW SUPPLIED

Item #	Description
8400ZA	QUINTIP® Individual Skin Test Device (sterile), case of 1800 units. Manufactured in England.
8406ZA	ComforTen® Covered 30-Hole Tray (Includes Spacers and Labels). Manufactured in the U.S.A.
8407ZA	ComforTen® Covered 60-Hole Tray (Includes Spacers and Labels). Manufactured in the U.S.A.
8924ZA	HollisterStier Skin Test Reservoirs (sterile), 30 per pouch. Manufactured in England.

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Jubilant HollisterStier LLC